

Biovek®



**Dyed (violet)
or undyed (beige)
braided multifilament
coated polyglycolide
(PGA)**

Synthetic
Absorbable



Suggested procedures

- Plastic
- Cosmetic
- Subcuticular
- Dental
- General / Intestinal
- Paediatric
- Urological
- Ophthalmic

Features















- Woven fine braid for smoother passage and minimal tissue trauma
- Optimum wound support
- Superior knot run down and knotting capabilities
- Coated with polycaprolactone and calcium stearate for a uniform absorption profile
- Absorption through hydrolysis provides predictable absorption irrespective of patient health

Tensile Strength	70%	35%
Days	14	21
Complete Absorption	60 - 90 days	

- Available in:
 - a range of needles including heavy needles
 - a range of micro needles (6mm to 10mm)
 - sizes 8/0 to USP 2
 - sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life

Commonly used Biovek® product codes and specifications.














Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	Colour	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
3/8 Circle Reverse Cutting									
10mm 	45cm	Violet	NAB601	NAB501					
10mm 	45cm	Violet	NADB601	NADB501					
12mm 	45cm	Violet	NAB602	NAB502					
13mm 	45cm	Violet	NAB6040	NAB5040	NAB4040				
		Undyed	NAB6040C	NAB5040C	NAB4040C				
35mm 	45cm	Violet					NAB2031	NAB031	NAB131
	76cm						NAB12031	NAB1031	NAB1131
40mm 	45cm	Violet					NAB2013	NAB013	NAB113
	76cm						NAB12013	NAB1013	NAB1113
45mm 	76cm	Violet					NAB12010	NAB1010	NAB1110
3/8 Circle Reverse Premium Cutting									
11mm 	45cm	Violet	NAB602D	NAB502D					
		Undyed	NAB602DC	NAB502DC					
13mm 	45cm	Violet	NAB6049	NAB5049					
		Undyed	NAB6049C	NAB5049C					
16mm 	45cm	Violet	NAB604	NAB504	NAB404	NAB304			
		Undyed		NAB504C	NAB404C	NAB304C			
	76cm	Violet		NAB1504	NAB1404	NAB1304			
		Undyed		NAB1504C	NAB1404C	NAB1304C			
19mm 	45cm	Violet		NAB505	NAB405	NAB305			
		Undyed		NAB505C	NAB405C	NAB305C			
	76cm	Violet		NAB1505	NAB1405	NAB1305			
		Undyed		NAB1505C	NAB1405C	NAB1305C			
24mm 	45cm	Violet			NAB406	NAB306	NAB206		
		Undyed			NAB406C	NAB306C	NAB206C		
	76cm	Violet			NAB1406	NAB1306	NAB1206		
		Undyed			NAB1406C	NAB1306C	NAB1206C		
26mm 	45cm	Violet				NAB3047	NAB2047		
		Undyed				NAB3047C	NAB2047C		
	76cm	Violet				NAB13047	NAB12047		
		Undyed				NAB13047C	NAB12047C		
30mm 	76cm	Violet					NAB1209		
		Undyed					NAB1209C		

* Needle images and sizes are approximations only.

Commonly used Biovek[®] product codes and specifications.










Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	Colour	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
1/2 Circle Reverse Cutting									
25mm 	76cm	Violet				NAB13034			
30mm 	45cm	Violet				NAB3038	NAB2038		
	76cm					NAB13038	NAB12038		
35mm 	76cm	Violet				NAB13017			
40mm 	76cm	Violet					NAB12012	NAB1012	NAB1112
	90cm						NAB2012S	NAB012S	NAB112S
48mm 	76cm	Violet						NAB1045	
	90cm							NAB045S	NAB145S
1/2 Circle Reverse Cutting (Heavy)									
35mm 	45cm	Violet					NAB2017E5	NAB017E5	
3/8 Circle Conventional Cutting									
13mm 	45cm	Violet	NAB6093	NAB5093	NAB4093				
16mm 	45cm	Undyed	NAB6037C	NAB5037C	NAB4037C				
24mm 	76cm	Undyed				NAB13023C	NAB12023C		
1/2 Circle Conventional Cutting									
22mm 	45cm	Violet				NAB3018			
48mm 	45cm	Undyed							NAB150C
Straight Premium Cutting									
<input checked="" type="checkbox"/> 	76cm	Violet					NAB2030		
60mm	90cm							NAB1030	
3/8 Circle Finline[®] (Conventional Premium) Cutting									
16mm 	45cm	Undyed	NAB6025C	NAB5025C	NAB4025C				
	76cm			NAB15025C					

* Needle images and sizes are approximations only.

Commonly used Biovek® product codes and specifications.








Contact us for any other suture specification.

Needle*	USP Metric	Colour	Thread								
			6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	
3/8 Circle Round Bodied Taper											
10mm 	76cm	Violet		NAB5052							
18mm 	76cm	Violet		NAB5066	NAB4066	NAB3066					
25mm 	76cm	Violet			NAB4067	NAB3067	NAB2067	NAB067			
		Undyed				NAB3067C	NAB2067C				
40mm 	76cm	Violet							NAB175	NAB275	
	90cm								NAB1175	NAB1275	
1/2 Circle Round Bodied Taper											
10mm 	76cm	Violet	NAB6051								
13mm 	76cm	Violet	NAB6055	NAB5055							
		Undyed	NAB6055C	NAB5055C							
16mm 	76cm	Violet		NAB5060	NAB4060	NAB3060	NAB2060				
18mm 	76cm	Violet		NAB5064	NAB4064	NAB3064	NAB2064				
20mm 	76cm	Violet		NAB5065	NAB4065	NAB3065	NAB2065				
22mm	76cm	Violet		NAB50#23	NAB40#23	NAB30#23	NAB20#23				
25mm	76cm	Violet			NAB4068	NAB3068	NAB2068	NAB068			
		Undyed				NAB3068C					
30mm	76cm	Violet			NAB4069	NAB3069	NAB2069	NAB069	NAB169	NAB269	
	90cm	Undyed					NAB2069C				
35mm	76cm	Violet				NAB3070	NAB2070	NAB070	NAB170		
	90cm	Undyed					NAB3070C	NAB2070C			
40mm	76cm	Violet						NAB2097	NAB097	NAB197	NAB297
	90cm							NAB12097	NAB1097	NAB1197	NAB1297
45mm	76cm	Violet						NAB20#15	NAB0#15	NAB1#15	NAB2#15
48mm	76cm	Violet						NAB2098	NAB098	NAB198	NAB298
	90cm							NAB12098	NAB1098	NAB1198	NAB1298

* Needle images and sizes are approximations only.

Commonly used Biovek® product codes and specifications.

















Contact us for any other suture specification.

Needle*	USP Metric	Colour	Thread				
			3/0 2	2/0 3	0 3.5	1 4	2 5
1/2 Circle Round Bodied Taper (Heavy)							
 35mm	76cm	Violet	NAB3070E5	NAB2070E5	NAB070E5	NAB170E5	
	90cm		NAB13070E5	NAB12070E5	NAB1070E5	NAB1170E5	
 40mm	76cm	Violet		NAB2097E5	NAB097E5	NAB197E5	NAB297E5
	90cm			NAB12097E5	NAB1097E5	NAB1197E5	NAB1297E5
 48mm	76cm	Violet		NAB2098E5	NAB098E5	NAB198E5	NAB298E5
	90cm			NAB12098E5	NAB1098E5	NAB1198E5	NAB1298E5
1/2 Circle Round Bodied Blunt Taper							
 40mm	76cm	Violet			NAB094	NAB194	
 48mm	90cm	Violet				NAB1196	
 64mm	90cm	Violet				NAB1182	
1/2 Circle Round Bodied Blunt Taper (Heavy)							
 48mm	90cm	Violet				NAB1196E5	
1/2 Circle Round Bodied A-CUTE® (Cutting Tip) Taper							
 18mm	76cm	Violet	NAB3035	NAB2035			
		Undyed	NAB3035C	NAB2035C			
 25mm	76cm	Violet	NAB3032	NAB2032			
 30mm	76cm	Violet		NAB2039	NAB039		

* Needle images and sizes are approximations only.

Commonly used Biovek® product codes and specifications.















Contact us for any other suture specification.

Needle*	USP Metric	Colour	Thread						
			4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	
1/2 Circle Round Bodied A-CUTE® (Cutting Tip) Taper									
 35mm	76cm	Violet			NAB2042	NAB042			
	90cm	Violet			NAB12042	NAB1042	NAB1142		
		Undyed			NAB12042C	NAB1042C	NAB1142C		
 40mm	90cm	Violet			NAB12079	NAB1079	NAB1179	NAB1279	
 45mm	76cm	Violet			NAB2044	NAB044	NAB144		
	90cm				NAB12044	NAB1044	NAB1144		
1/2 Circle Round Bodied A-CUTE® (Cutting Tip) Taper (Heavy)									
 35mm	76cm	Violet			NAB2042E5	NAB042E5	NAB142E5		
	90cm	Undyed					NAB1142E5C		
Sterile Pre-Cut Lengths - Non-Needle									
N/A		1x100cm	Violet	ABA4040	ABA3040	ABA2040	ABA040	ABA140	
N/A		1x150cm	Violet	ABA4060	ABA3060	ABA2060	ABA060	ABA160	ABA260
N/A		1x180cm	Violet	ABA4072	ABA3072	ABA2072	ABA072	ABA172	
N/A		1x250cm	Violet	ABA40100	ABA30100	ABA20100			
N/A		2x76cm	Violet	ABB4030	ABB3030	ABB2030	ABB030		
N/A		5x60cm	Violet	ABE4024	ABE3024	ABE2024	ABE024	ABE124	ABE224
N/A		5x76cm	Violet	ABE4030	ABE3030	ABE2030	ABE030	ABE130	
N/A		6x45cm	Violet	ABF4018	ABF3018	ABF2018	ABF018	ABF118	
N/A		10x60cm	Violet	ABJ4024	ABJ3024	ABJ2024	ABJ024	ABJ124	
N/A		12x45cm	Violet	ABL4018	ABL3018	ABL2018	ABL018	ABL118	
N/A		12x60cm	Violet	ABL4024	ABL3024	ABL2024	ABL024	ABL124	
N/A		13x60cm	Violet	ABM4024	ABM3024	ABM2024	ABM024	ABM124	ABM224

* Needle images and sizes are approximations only.

Commonly used Biovek® product codes and specifications.

Contact us for any other suture specification.

Micro Suture Range		Thread					
Needle*	USP Metric	Colour	8/0 0.4	7/0 0.5	6/0 0.7	5/0 1	
1/4 Circle Lancet Spatula							
8.0mm (wire size 0.30mm)  	45cm	Violet	ADB8F8020	ADB7F8020	ADB6F8020	ADB5F8020	
8.0mm (wire size 0.35mm)  	45cm	Violet	ADB8F8021	ADB7F8021	ABD6F8021	ADB5F8021	
8.5mm (wire size 0.30mm)  	45cm	Violet	ADB8F8520	ADB7F8520	ADB6F8520	ADB5F8520	
3/8 Circle Lancet Spatula							
6.0mm (wire size 0.20mm)  	45cm	Violet	ADB8F6011	ADB7F6011	ADB6F6011		
6.4mm (wire size 0.20mm)  	45cm	Violet	ADB8F6411	ADB7F6411	ADB6F6411	ADB5F6411	
8.0mm (wire size 0.20mm)  	45cm	Violet		ADB7F8011	ADB6F8011	ADB5F8011	
1/2 Circle Lancet Spatula							
8.0mm (wire size 0.20mm)  	45cm	Violet		ADB7F8013	ADB6F8013	ADB5F8013	
8.0mm (wire size 0.25mm)  	45cm	Violet		ADB7F8023	ADB6F8023	ADB5F8023	
3/8 Circle Reverse Cutting							
6.4mm (wire size 0.28mm)  	45cm	Violet		ADB7F6414	ADB6F6414		

* Needle images and sizes are approximations only.

Proudly Australian owned and operated.

Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



[dynek.com](https://www.dynek.com)

AU: 1800 088 172

NZ: 0800 459 649

INT: +61 8 8268 2033

sales@dynek.com



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dynek sutures.
- Biovek® is included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.

QA:BE063/4 06/2022