

# Nylene<sup>®</sup>



## Blue monofilament polyamide (Nylon)

Synthetic  
Non-absorbable



### Suggested procedures












- General
- Dermatology
- Plastic
- Cosmetic
- Accident / Emergency
- Subcuticular
- Skin cancer
- General / Intestinal

### Features

- Physiologically inert with minimal tissue reaction
- Pliable, smooth, easy to handle, with excellent knot holding ability
- Extruded through a diamond die into a single uniform monofilament strand
- Extremely strong with high tensile strength
- Bright blue colour provides increased visibility
- Also available in black and clear (undyed) for some suture specifications
- Also available in sterile pre-cut lengths
- Sutures made with 300 series stainless steel
- Up to five-year shelf life

# Commonly used Nylene® product codes and specifications.










Contact us for any other suture specification.

Needle*		Thread									
		USP Metric	Colour	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
<b>3/8 Circle Reverse Cutting</b>											
10mm		45cm	Blue	N601	N501						
12mm		45cm	Blue	N602	N502	N402					
13mm		45cm	Blue	N6040	N5040						
15mm		45cm	Blue			N404E10E20					
20mm		45cm	Blue		N5046						
28mm		45cm	Blue					N207	N07		
		76cm	Blue					N1207			
35mm		45cm	Blue			N4031	N3031	N2031			
			Clear			N4031C					
		76cm	Blue			N14031	N13031	N12031			
40mm		45cm	Blue				N3013	N2013	N013	N113	N213
		76cm					N13013	N12013	N1013	N1113	N1213
		90cm					N3013S	N2013S			
45mm		45cm	Blue					N2010	N010	N110	
		76cm						N12010	N1010	N1110	
77mm		45cm	Blue							N141	N241
77mm		76cm	Blue							DN1141	DN1241

\* Needle images and sizes are approximations only.

# Commonly used Nylene® product codes and specifications.













Contact us for any other suture specification.

Needle*		Thread								
		USP Metric	Colour	7/0 0.5	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
<b>3/8 Circle Reverse Premium Cutting</b>										
10.5mm	<input checked="" type="checkbox"/> 	45cm	Blue	N702DE10	N602DE10					
11mm	<input checked="" type="checkbox"/> 	45cm	Blue	N702D	N602D	N502D				
13mm	<input checked="" type="checkbox"/> 	45cm	Blue		N6049	N5049	N4049			
		76cm	Black		N6049B	N5049B				
16mm	<input checked="" type="checkbox"/> 	45cm	Blue		N604	N504	N404	N304		
		76cm	Black		N604B	N504B				
		90cm	Blue		N1604	N1504	N1404	N1304		
19mm	<input checked="" type="checkbox"/> 	45cm	Blue		N605	N505	N405	N305	N205	
			Clear				N405C	N305C		
		76cm	Blue		N1605	N1505	N1405	N1305	N1205	
		90cm	Clear				N1405C	N1305C		
20mm	<input checked="" type="checkbox"/> 	45cm	Blue			N503	N403	N303		
			Clear							
24mm	<input checked="" type="checkbox"/> 	45cm	Blue		N606	N506	N406	N306	N206	
			Black				N406B	N306B		
		76cm	Blue		N1606	N1506	N1406	N1306	N1206	
		90cm	Black					N1306B	N1206B	
26mm	<input checked="" type="checkbox"/> 	45cm	Blue					N3047	N2047	
		76cm	Blue					N13047	N12047	
30mm	<input checked="" type="checkbox"/> 	45cm	Blue					N309	N209	N09
		76cm	Blue					N1309	N1209	N109
			Black					N1309B	N1209B	

\* Needle images and sizes are approximations only.

# Commonly used Nylene® product codes and specifications.











Contact us for any other suture specification.

Needle*		Thread								
		USP Metric	Colour	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
<b>1/2 Circle Reverse Cutting</b>										
16mm		45cm	Blue	N5015	N4015	N3015				
18mm		45cm	Blue	N5014	N4014	N3014				
		76cm		N15014	N14014	N13014				
25mm		45cm	Blue		N4034	N3034	N2034			N234
		76cm			N14034	N13034	N12034			N1234
30mm		45cm	Blue			N3038	N2038			
		76cm			N14038	N13038				
35mm		45cm	Blue			N3017	N2017	N017	N117	N217
		76cm				N13017	N12017	N1017	N1117	N1217
40mm		45cm	Blue				N2012	N012	N112	N212
		76cm					N12012	N1012	N1112	N1212
48mm		76cm	Blue							N1245
		90cm						N045S	N145S	N245S
<b>1/2 Circle Reverse Premium Cutting</b>										
15mm		45cm	Blue	N5016	N4016					
		76cm	Black		N4016B					
		76cm	Blue		N14016					
<b>Straight Premium Cutting</b>										
60mm		100cm	Blue			N3030	N2030			
<b>3/8 Circle Conventional Cutting</b>										
16mm		45cm	Black	N5037B						
40mm		45cm	Blue				N2043			
90mm		50cm	Black						N111BL50	
		76cm	Blue				N2011	N011	N111	N211
		100cm	Blue				N12011	N1011	N1111	N1211

\* Needle images and sizes are approximations only.

# Commonly used Nylene® product codes and specifications.

Contact us for any other suture specification.

		Thread								
Needle*	USP Metric	Colour	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
<b>3/8 Circle Finline® (Conventional Premium) Cutting</b>										
16mm 	45cm	Blue	N6025	N5025	N4025	N3025				
		Clear	N6025C	N5025C	N4025C	N3025C				
20mm 	45cm	Blue		N5026	N4026	N3026				
		Clear		N5026C	N4026C	N3026C				
25mm 	45cm	Blue		N5027	N4027	N3027				
		Clear		N5027C	N4027C	N3027C				
<b>3/8 Circle Round Bodied Taper</b>										
13mm 	76cm	Blue		N5053						
16mm 	76cm	Blue			N4054	N3054				
<b>1/2 Circle Round Bodied Taper</b>										
18mm 	76cm	Blue		N5064	N4064					
25mm 	76cm	Blue			N4068	N3068	N2068	N068		
	100cm							N1068	N1168	
30mm 	76cm	Blue				N3069	N2069	N069	N169	
	100cm							N1069	N1169	
35mm 	76cm	Blue					N2070	N070	N170	
40mm 	100cm	Blue					N2097	N097	N197	
		Black							N197B	
48mm 	76cm	Blue							N198	
	100cm								N1198	N1298

\* Needle images and sizes are approximations only.

# Commonly used Nylene® product codes and specifications.

Contact us for any other suture specification.

Needle*	Thread					
	USP Metric	Colour	2/0 3	0 3.5	1 4	2 5
<b>1/2 Circle Round Bodied Taper (Heavy)</b>						
35mm 	76cm	Blue	N2070E5	N070E5	N170E5	
40mm 	100cm	Blue	N2097E5	N097E5	N197E5	N297E5
48mm 	76cm	Blue	N2098E5	N098E5	N198E5	N298E5
	100cm				N1198E5	N1298E5
<b>1/2 Circle Round Bodied Blunt Point</b>						
35mm 	100cm	Blue			N11#13	
48mm 	150cm	Blue			N196L150	
<b>1/2 Circle Round Bodied Taper (LOOP)</b>						
35mm 	150cm	Blue			N170LOOP	
40mm 	150cm	Blue		N097LOOP	N197LOOP	N297LOOP
48mm 	140cm	Blue			N198LOOP140	
	150cm	Blue			N198LOOP	N298LOOP
	200cm	Blue			N198LOOP200	

\* Needle images and sizes are approximations only.

# Commonly used Nylene® product codes and specifications.

Contact us for any other suture specification.

Needle*		Thread									
		USP Metric	Colour	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
Sterile Pre-Cut Lengths - Non-Needle											
N/A		1x76cm	Blue	NA6030	NA5030	NA4030	NA3030	NA2030	NA030	NA130	NA230
N/A		1x100cm	Blue	NA6040	NA5040	NA4040	NA3040	NA2040	NA040	NA140	NA240
N/A		1x180cm	Blue	NA6072	NA5072	NA4072	NA3072	NA2072	NA072	NA172	NA272
N/A		2x76cm	Blue	NB6030	NB5030	NB4030	NB3030	NB2030	NB030	NB130	NB230
N/A		2x100cm	Blue	NB6040	NB5040	NB4040	NB3040	NB2040	NB040	NB140	NB240
N/A		2x180cm	Blue	NB6072	NB5072	NB4072	NB3072	NB2072	NB072	NB172	NB272
N/A		3x76cm	Blue	NC6030	NC5030	NC4030	NC3030	NC2030	NC030	NC130	NC230
N/A		3x100cm	Blue	NC6040	NC5040	NC4040	NC3040	NC2040	NC040	NC140	NC240
N/A		3x180cm	Blue	NC6072	NC5072	NC4072	NC3072	NC2072	NC072	NC172	NC272

\* Needle images and sizes are approximations only.

# Proudly Australian owned and operated.

Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



[dynek.com](https://dynek.com)

**AU: 1800 088 172**

**NZ: 0800 459 649**

**INT: +61 8 8268 2033**

[sales@dynek.com](mailto:sales@dynek.com)



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dynek sutures.
- Nylene® is included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.