

# Surgical Catgut



## Monofilament plain and chromic surgical catgut

Natural  
Absorbable

### Suggested procedures

- Dental
- Subcuticular
- Plastics
- General / Intestinal
- Obstetrics & gynaecology
- Urological
- Oral

### Features

- A specialised spinning and grinding process creates a suture strand of uniform diameter
- The suture strand is further polished to achieve maximum smoothness for reliability and strength
- Uniform diameter of plain and chromic gut allows for an even absorption rate
- Chromic catgut is chromicised before it is spun to maintain the even absorption rate
- Optimal tensile strength and in vivo performance
- Manufactured from regulatory approved intestinal serosa of cattle
- Absorption through enzyme digestion:


















<b>Plain</b>	7-10 days
<b>Chromic</b>	20-40 days

- Available in sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life



# Commonly used Chromic Catgut product codes and specifications.

















Contact us for any other suture specification.

		Thread							
Needle*	USP Metric	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	
<b>3/8 Circle Reverse Cutting</b>									
10.5mm		50cm	DCG601E10						
12mm		50cm	CG602	CG502	CG402				
12mm		50cm	DCG602						
13mm		30cm		CG5040L30	CG4040L30				
		50cm		CG5040	CG4040				
<b>3/8 Circle Reverse Premium Cutting</b>									
11mm		50cm	CG602D	CG502D					
13mm		50cm	CG6049	CG5049	CG4049				
16mm		50cm		CG504	CG404	CG304	CG204		
		76cm		CG1504	CG1404	CG1304	CG1204		
19mm		50cm		CG505	CG405	CG305	CG205		
		76cm		CG1505	CG1405	CG1305	CG1205		
24mm		50cm			CG406	CG306	CG206		
		76cm			CG1406	CG1306	CG1206		
26mm		50cm			CG4047	CG3047	CG2047		
30mm		50cm				CG309	CG209		
		76cm				CG1309	CG1209		
<b>1/2 Circle Reverse Cutting</b>									
16mm		50cm		CG5015	CG4015				
25mm		76cm			CG14034	CG13034	CG12034	CG1134	
35mm		76cm				CG3017	CG2017	CG017	CG117
		100cm				CG13017	CG12017	CG1017	CG1117
40mm		76cm					CG1012	CG1112	
48mm		50cm					CG045	CG145	
<b>3/8 Circle Finline® (Conventional Premium) Cutting</b>									
16mm		50cm			CG4025				

\* Needle images and sizes are approximations only.

# Commonly used Chromic Catgut product codes and specifications.











Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
<b>3/8 Circle Conventional Cutting</b>									
13mm 	50cm	CG6093	CG5093	CG4093					
16mm 	50cm	CG6037	CG5037	CG4037	CG3037				
<b>3/8 Circle Round Bodied Taper</b>									
16mm 	76cm			CG4054	CG3054	CG2054			
25mm 	76cm				CG3067	CG2067			
40mm 	76cm					CG2075			
<b>1/2 Circle Round Bodied Taper</b>									
16mm 	76cm			CG4060	CG3060				
18mm 	76cm		CG5064	CG4064	CG3064	CG2064			
20mm 	76cm		CG5065	CG4065	CG3065	CG2065			
22mm 	76cm			CG40#23	CG30#23				
25mm 	76cm		CG5068	CG4068	CG3068	CG2068	CG068	CG168	
25mm 	76cm				DCG3068	DCG2068			
30mm 	76cm				CG3069	CG2069	CG069	CG169	
35mm 	76cm				CG3070	CG2070	CG070	CG170	CG270
	100cm				CG13070	CG12070	CG1070	CG1170	
40mm 	76cm					CG2097	CG097	CG197	CG297
	100cm					CG12097	CG1097	CG1197	CG1297
48mm 	76cm						CG098	CG198	CG298
	100cm						CG1098	CG1198	CG1298
65mm 	76cm					CG2085		CG185	

\* Needle images and sizes are approximations only.

# Commonly used Chromic Catgut product codes and specifications.

















Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6	3 7	
<b>1/2 Circle Round Bodied Taper</b>									
 76mm	76cm				CG2086	CG086			
<b>1/2 Circle Round Bodied Taper &amp; Reverse Cutting</b>									
 35mm	150cm				DCG2070E1				
<b>1/2 Circle Round Bodied Taper (Heavy)</b>									
 35mm	76cm				CG2070E5	CG070E5	CG170E5		
 40mm	76cm				CG2097E5	CG097E5	CG197E5	CG297E5	
	100cm				CG12097E5	CG1097E5	CG1197E5	CG1297E5	CG1397E5
 48mm	76cm					CG098E5	CG198E5	CG298E5	
	100cm					CG1098E5	CG1198E5	CG1298E5	
<b>1/2 Circle Round Bodied A-CUTE® (Cutting Tip) Taper</b>									
 35mm	76cm				CG2042	CG042			
 45mm	76cm				CG2044	CG044	CG144	CG244	
<b>3/8 Circle Round Bodied Blunt Point</b>									
 64mm	76cm						CG182		
<b>Sterile Pre-Cut Lengths - Non-Needle</b>									
N/A		150cm	CGA4060	CGA3060	CGA2060	CGA060	CGA160	CGA260	CGA360
N/A		250cm	CGA40100	CGA30100	CGA20100	CGA0100	CGA1100	CGA2100	CGA3100

\* Needle images and sizes are approximations only.

# Commonly used Plain Catgut product codes and specifications.













Contact us for any other suture specification.

Needle*	Thread							
	USP Metric	6/0 1	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5
<b>3/8 Circle Reverse Cutting</b>								
12mm 	50cm	PG602	PG502	PG402				
<b>3/8 Circle Reverse Premium Cutting</b>								
10.5mm 	50cm	PG602DE10						
 10.5mm 	50cm	DPG602DE10						
11mm 	50cm	PG602D	PG502D					
 11mm 	50cm	DPG602D						
13mm 	50cm	PG6049	PG5049	PG4049				
16mm 	50cm		PG504	PG404	PG304			
	76cm		PG1504	PG1404	PG1304			
19mm 	50cm		PG505	PG405	PG305	PG205		
	76cm		PG1505	PG1405	PG1305	PG1205		
24mm 	50cm			PG406	PG306	PG206		
	76cm			PG1406	PG1306	PG1206		
26mm 	50cm				PG3047	PG2047		
	76cm				PG13047	PG12047		
30mm 	50cm				PG309	PG209		
	76cm				PG1309	PG1209		
<b>1/2 Circle Reverse Cutting</b>								
16mm 	76cm				PG13015			
25mm 	50cm				PG3034			
40mm 	76cm						PG1012	PG1112

\* Needle images and sizes are approximations only.

# Commonly used Plain Catgut product codes and specifications.








Contact us for any other suture specification.

Needle*	Thread						
	USP Metric	5/0 1.5	4/0 2	3/0 3	2/0 3.5	0 4	1 5
<b>3/8 Circle Conventional Cutting</b>							
16mm 	76cm		PG14037				
<b>Straight Premium Cutting</b>							
<input checked="" type="checkbox"/>  60mm	100cm			PG3030	PG2030		
<b>3/8 Circle Round Bodied Taper</b>							
16mm 	76cm		PG4054				
20mm 	76cm			PG3057			
<b>1/2 Circle Round Bodied Taper</b>							
18mm 	76cm		PG4064				
20mm 	76cm	PG5065	PG4065	PG3065	PG2065		
25mm 	76cm		PG4068	PG3068	PG2068	PG068	PG168
30mm 	76cm			PG3069	PG2069	PG069	PG169
35mm 	76cm			PG3070	PG2070	PG070	PG170
40mm 	76cm				PG2097	PG097	PG197
48mm 	76cm				PG2098	PG098	PG198
76mm 	76cm				PG2086	PG086	

\* Needle images and sizes are approximations only.

# Commonly used Plain Catgut product codes and specifications.

Contact us for any other suture specification.

Needle*		Thread						
		USP Metric	4/0 2	3/0 3	2/0 3.5	0 4	1 5	2 6
<b>1/2 Circle Round Bodied Taper (Heavy)</b>								
35mm		76cm					PG070E5	PG170E5
40mm		76cm			PG2097E5	PG097E5	PG197E5	PG297E5
48mm		76cm			PG2098E5	PG098E5	PG198E5	PG298E5
		100cm			PG12098E5	PG1098E5	PG1198E5	PG1298E5
<b>1/2 Circle Round Bodied A-CUTE® (Cutting Tip) Taper</b>								
25mm		50cm		PG3032Z	PG2032Z			
		76cm		PG3032	PG2032			
35mm		76cm			PG2042			
<b>Sterile Pre-Cut Lengths - Non-Needle</b>								
N/A		150cm	PGA4060	PGA3060	PGA2060	PGA060	PGA160	PGA260
N/A		250cm	PGA40100	PGA30100	PGA20100	PGA0100	PGA1100	PGA2100

\* Needle images and sizes are approximations only.

# Proudly Australian owned and operated.

Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.




**dynek.com**

**AU: 1800 088 172**

**NZ: 0800 459 649**

**INT: +61 8 8268 2033**

**sales@dynek.com**



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dynek sutures.
- Dynek surgical catgut is included on the ARTG (AU) and Medsafe (NZ).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.

QA:BE080/1 06/2022