

# Dyflex<sup>®</sup> Teflex<sup>®</sup>



**Green or white  
braided multifilament  
coated polyester**

Synthetic  
Non-absorbable



## **Suggested procedures**
















- Subcuticular
- Orthopaedics
- Plastic
- General / Intestinal
- Obstetrics

## **Features**

- Excellent handling and tying properties
- Very low memory prevents knotting and tangling during suture removal
- Silicone treated to minimise tissue reaction
- Minimal stretch maintaining tensile strength
- Uniform diameter throughout
- Dyflex<sup>®</sup> has been treated to increase softness and smoothness
- Teflex<sup>®</sup> has been designed to allow its use as a subcuticular suture
- Available in:
  - single-armed and double-armed specifications
  - various multi-packs
  - a range of needles including heavy needles
  - sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life

# Commonly used Dyflex<sup>®</sup> product codes and specifications.

Contact us for any other suture specification.
















Needle*		Thread								
		USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
<b>3/8 Circle Reverse Cutting</b>										
12mm		45cm	P602S	P502S						
35mm		76cm					P12031S	P1031S	P1131S	
40mm		76cm					P12013S	P1013S	P1113S	
<b>3/8 Circle Reverse Premium Cutting</b>										
13mm		45cm		P5049S						
		76cm		P15049S						
16mm		45cm		P504S	P404S					
19mm		45cm		P505S	P405S	P305S				
24mm		45cm			P406S	P306S	P206S			
26mm		45cm				P3047S	P2047S			
30mm		76cm					P1209S	P109S		
<b>1/2 Circle Reverse Cutting</b>										
20mm		45cm								P287S
22mm		76cm							P11#88S	P12#88S
25mm		45cm							P134S	
		76cm							P1134S	P1234S
35mm		76cm						P1017S	P1117S	
40mm		76cm						P1012S	P1112S	P1212S
<b>1/2 Circle Reverse Cutting (Loop)</b>										
60mm		150cm								P290SLOOP

\* Needle images and sizes are approximations only.

# Commonly used Dyflex<sup>®</sup> product codes and specifications.

Contact us for any other suture specification.














Needle*	Thread						
	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5
<b>Straight Premium Cutting</b>							
 60mm	100cm					P3030S	P2030S
<b>3/8 Circle Round Bodied Taper</b>							
 13mm	76cm	P6053S	P5053S	P4053S			
 13mm	76cm	DP6053S	DP5053S				
 16mm	76cm		P5054S	P4054S	P3054S		
 16mm	76cm		DP5054S	DP4054S			
<b>1/2 Circle Round Bodied Taper</b>							
 13mm	76cm	P6055S	P5055S	P4055S			
 13mm	76cm	DP6055S	DP5055S	DP4055S			
 16mm	76cm		P5060S	P4060S	P3060S		
 16mm	76cm		DP5060S	DP4060S	DP3060S	DP2060S	
 18mm	76cm		P5064S	P4064S	P3064S	P2064S	
 18mm	76cm	DP6064S	DP5064S	DP4064S	DP3064S	DP2064S	
	90cm			DP14064S			
 20mm	76cm		P5065S	P4065S	P3065S	P2065S	
 20mm	76cm		DP5065S	DP4065S	DP3065S	DP2065S	
 25mm	76cm			P4068S	P3068S	P2068S	P068S
	90cm					P12068S	P1068S
 25mm	76cm			DP4068S	DP3068S	DP2068S	DP068S

\* Needle images and sizes are approximations only.

# Commonly used Dyflex<sup>®</sup> product codes and specifications.
















Contact us for any other suture specification.

Needle*	Thread									
	USP Metric	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	3+4 6	5 7
<b>1/2 Circle Round Bodied Taper</b>										
 30mm	76cm			P3069S	P2069S	P069S	P169S			
 30mm	76cm			DP3069S	DP2069S	DP069S				
 35mm	76cm			P3070S	P2070S	P070S	P170S	P270S		P570S
 35mm	76cm 90cm			DP3070S	DP2070S	DP070S				
 40mm	76cm								P297S	
 40mm	76cm						DP097S	DP197S	DP297S	
 48mm	76cm							P198S		
<b>1/2 Circle Round Bodied Taper (Heavy)</b>										
 48mm	76cm						P198SE5	P298SE5	P398SE5	P598SE5
<b>1/2 Circle Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper</b>										
 18mm	76cm	P5035S	P4035S	P3035S	P2035S					
 18mm	76cm	DP5035S	DP4035S	DP3035S	DP2035S					
 25mm	76cm		P4032S	P3032S	P2032S					

\* Needle images and sizes are approximations only.

# Commonly used Dyflex<sup>®</sup> product codes and specifications.

Contact us for any other suture specification.

		Thread							
		USP Metric	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	5 7
<b>1/2 Circle Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper</b>									
25mm		76cm	DP4032S	DP3032S	DP2032S				
		90cm	DP14032S	DP13032S	DP12032S	DP1032S			
40mm		76cm				P079S			
55mm		76cm						P583S	
<b>1/2 Circle Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper (Heavy)</b>									
45mm		90cm						P244SMDE5	
<b>4/5 Curve Round Bodied Taper J Needle</b>									
16mm		76cm					P1#38S		
<b>4/5 Curve Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper J Needle</b>									
30mm		76cm					P1#42S		
<b>Sterile Pre-Cut Lengths - Non-Needle</b>									
N/A		1x76cm						PA530S	
		1x100cm					PA140S	PA240S	
N/A		1x150cm		PA3060S	PA2060S				
N/A		1x180cm			PA2072S	PA072S	PA172S	PA272S	
N/A		1x200cm		PA3080S	PA2080S	PA080S	PA180S		
N/A		1x250cm		PA30100S	PA20100S	PA0100S	PA1100S		
N/A		1x300cm	PA40120S	PA30120S	PA20120S	PA0120S	PA1120S		
N/A		3x45cm						PC218S	
N/A		6x45cm				PF018S			
N/A		6x60cm						PF224S	

\* Needle images and sizes are approximations only.

# Commonly used Teflex<sup>®</sup> product codes and specifications.

Contact us for any other suture specification.










Needle*		Thread				
		USP Metric	5/0 1	4/0 1.5	3/0 2	1 4
<b>3/8 Circle Reverse Premium Cutting</b>						
13mm		45cm		T4049		
		76cm		T14049		
19mm		45cm	T505	T405	T305	
<b>1/2 Circle Reverse Cutting</b>						
25mm		76cm				T1134 T1234
40mm		45cm				T112
<b>3/8 Circle Finesline<sup>®</sup> (Conventional Premium) Cutting</b>						
16mm		45cm	T5025	T4025	T3025	
16mm		45cm	DT5025	DT4025	DT3025	

\* Needle images and sizes are approximations only.

# Commonly used Dyflex<sup>®</sup> and Teflex<sup>®</sup> Multi Pack codes and specifications

Contact us for any other suture specification.



Needle*	Thread							
	USP Metric	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5	5 7
<b>1/2 Circle Round Bodied Taper - Multi Pack - 5 Green &amp; 5 White pieces</b>								
18mm CV300 	90cm			DPT2064SMJ				
20mm CV300 	90cm			DPT2065SMJ				
25mm CV300 	90cm			DPT2068SMJ				
<b>1/2 Circle Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper - Multi Pack - 5 Green &amp; 5 White pieces</b>								
18mm CV300 	76cm			DPT2035SZMJ				
	90cm			DPT2035SMJ				
25mm CV300 	76cm			DPT2032SZMJ				
	90cm			DPT2032SMJ				
<b>1/2 Circle Round Bodied A-CUTE<sup>®</sup> (Cutting Tip) Taper - Multi Pack - 3 Green &amp; 3 White pieces</b>								
18mm CV300 	90cm			DPT2035SMF				
25mm CV300 	90cm			DPT2032SMF				

\* Needle images and sizes are approximations only.

# Proudly Australian owned and operated.

Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



[dynek.com](https://www.dynek.com)

**AU: 1800 088 172**

**NZ: 0800 459 649**

**INT: +61 8 8268 2033**

[sales@dynek.com](mailto:sales@dynek.com)



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dynek sutures.
- Dyflex® and Teflex® are included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.

QA: BE081/1 06/2022