


Dysilk[®]



Black braided multifilament silk

Natural
Non-absorbable



Suggested procedures





















- Dental
- Oral maxillo facial
- Ophthalmic
- Subcuticular
- General / Intestinal
- Plastic

Features

- Excellent handling and tying properties
- Very low memory prevents knotting and tangling during suture removal
- Specially treated braid allows ease of handling
- Braid made of a core of twisted silk around which a cover is plaited
- Uniform thread diameter and surface smoothness for good knotting characteristics
- Doesn't become hard or difficult to handle
- Ties down with a high degree of knot security
- Available in sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life

Commonly used Dysilk® product codes and specifications.






















Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
3/8 Circle Reverse Cutting									
10.5mm 	45cm	S601E10							
12mm 	45cm	S602	S502	S402					
	76cm		S1502						
13mm 	45cm	S6040		S4040					
28mm 	45cm				S307	S207	S07		
	76cm				S1307	S1207	S107		
35mm 	45cm					S2031			
	76cm				S13031	S12031			
40mm 	76cm					S2013	S1013	S1113	
45mm 	50cm					S2010	S010		
3/8 Circle Reverse Premium Cutting									
11mm 	45cm	DS602D							
13mm 	45cm		S5049						
16mm 	45cm	S604	S504	S404	S304	S204			
	76cm		S1504	S1404	S1304	S1204			
19mm 	45cm	S605	S505	S405	S305				
	76cm		S1505	S1405	S1305				
24mm 	45cm			S406	S306	S206			
	76cm			S1406	S1306	S1206			
26mm 	45cm				S3047	S2047			
30mm 	45cm				S309	S209	S09		S29
	76cm				S1309	S1209	S109		
1/2 Circle Reverse Cutting									
13mm 	45cm			S40#7					
16mm 	45cm			S4015					
18mm 	45cm			S4014	S3014				
25mm 	76cm					S12034			
30mm 	76cm				S13038	S12038	S1038		
40mm 	45cm						S012		
	76cm						S1012	S1112	

* Needle images and sizes are approximations only.

Commonly used Dysilk[®] product codes and specifications.

Contact us for any other suture specification.

		Thread							
Needle*	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
Straight Premium Cutting									
60mm 	100cm				S3030	S2030	S030		
3/8 Circle Conventional Cutting									
16mm 	45cm		S5037						
	76cm		S15037	S14037					
3/8 Circle Finline[®] (Conventional Premium) Cutting									
16mm 	45cm	S6025	S5025	S4025	S3025				
20mm 	45cm			S4026	S3026				
25mm 	45cm		S5027	S4027					
3/8 Circle Round Bodied Taper									
10mm 	76cm		S5052	S4052					
13mm 	76cm	S6053	S5053						
16mm 	76cm		S5054	S4054	S3054				
1/2 Circle Round Bodied Taper									
16mm 	76cm			S4060	S3060	S2060			
20mm 	76cm		S5065	S4065	S3065	S2065			
25mm 	76cm			S4068	S3068	S2068	S068	S168	
30mm 	76cm				S3069	S2069	S069	S169	
35mm 	76cm					S2070	S070	S170	
1/2 Circle Round Bodied A-CUTE[®] (Cutting Tip) Taper									
25mm 	76cm					S2032			
Sterile Pre-Cut Lengths - Non-Needle									
N/A 	1x100cm			SA4040	SA3040	SA2040	SA040	SA140	
N/A 	1x180cm			SA4072	SA3072	SA2072	SA072	SA172	
N/A 	5x60cm			SE4024	SE3024	SE2024	SE024	SE124	
N/A 	13x60cm			SM4024	SM3024	SM2024	SM024	SM124	
N/A 	6x45cm			SF4018	SF3018	SF2018	SF018	SF118	
N/A 	10x76cm			SJ4030	SJ3030	SJ2030	SJ030	SJ130	
N/A 	12x76cm					SL2030			

* Needle images and sizes are approximations only.

Proudly Australian owned and operated.

Dyneke is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dyneke has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



dyneke.com

AU: 1800 088 172

NZ: 0800 459 649

INT: +61 8 8268 2033

sales@dyneke.com



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dyneke sutures.
- Dysilk® is included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.